

MicroWorldOmics install tutorial

(Windows)

Download zip package:

For Chinese main land:

Download from links 1:

<https://caiyun.139.com/w/i/2oRhbaZRe2g7e>; password: 9pm8

Download from links 2:

<https://www.alipan.com/s/qZ7DFgDJx8i>

Download from links 3:

https://pan.baidu.com/s/1nGFqxZOJWSmlM_EqgHz7SA; password: smrz

For World users:


Download from links:

<https://drive.google.com/drive/folders/1DXQ9gOUdg4zVWxAa2w3wKoV7rT-4UJAc?usp=sharing>


Notice: The installation path must be in English and cannot contain Spaces.

Download zip package:


1. Decompress module package:

 MicroWorldOmics_modules.rar


for example, the decompression path is “D:\MicroWorldOmics\”




2. Double-click “Setup.exe”:



3. Select “I Agree” for the copyright:
















4. User-defined installation folder, click install and finish installation:



Note that the installation path should be consistent with the decompression path of the module package.

5. Check installation:

	aiohttp	2024/4/9 13:48	文件夹
	altgraph-0.17.3.dist-info	2024/4/9 13:48	文件夹
	bcrypt	2024/4/9 13:48	文件夹
	Bio	2024/4/9 13:48	文件夹
	blast-BLAST_VERSION+	2024/4/9 13:49	文件夹
	caffe2	2024/4/9 13:49	文件夹
	certifi	2024/4/9 13:49	文件夹
	coverage	2024/4/9 13:49	文件夹
	cryptography	2024/4/9 13:49	文件夹
	cryptography-40.0.2.dist-info	2024/4/9 13:49	文件夹
	cv2	2024/4/9 13:49	文件夹
	MicroWorldOmics.exe	2025/4/10 15:04	应用程序
	MicroWorldOmics	2025/6/26 10:26	Internet 快捷方式

The module package needs to be in the same directory as MicroWorldOmics.exe.

6. Open and launch “MicroWorldOmics” in your start menu:



MicroWorldOmics

应用



应用



Microsoft Store



Excel



PowerPoint



WPS Office

