

# MicroWorldOmics install tutorial

(Windows)

## Download zip package:

For Chinese main land, download from links:

<https://caiyun.139.com/m/i?2jQXoNqM9JV3r>; password: d3p9


For World users, download from links:

[https://drive.google.com/drive/folders/1dggfsdfXisMwy0\\_Xoo8LoGUqAmmOClaX?usp=sharing](https://drive.google.com/drive/folders/1dggfsdfXisMwy0_Xoo8LoGUqAmmOClaX?usp=sharing)


**Notice: The installation path must be in English and cannot contain Spaces.**

## Download zip package:


1. Decompress module package:




for example, the decompression path is “D:\MicroWorldOmics\”




2. Double-click “Setup.exe”:



3. Select “I Agree” for the copyright:
















4. User-defined installation folder, click install and finish installation:



Note that the installation path should be consistent with the decompression path of the module package.

5. Check installation:

	aiohttp	2024/4/9 13:48	文件夹
	altgraph-0.17.3.dist-info	2024/4/9 13:48	文件夹
	bcrypt	2024/4/9 13:48	文件夹
	Bio	2024/4/9 13:48	文件夹
	blast-BLAST_VERSION+	2024/4/9 13:49	文件夹
	caffe2	2024/4/9 13:49	文件夹
	certifi	2024/4/9 13:49	文件夹
	coverage	2024/4/9 13:49	文件夹
	cryptography	2024/4/9 13:49	文件夹
	cryptography-40.0.2.dist-info	2024/4/9 13:49	文件夹
	cv2	2024/4/9 13:49	文件夹
	MicroWorldOmics.exe	2025/4/10 15:04	应用程序
	MicroWorldOmics	2025/6/26 10:26	Internet 快捷方式

The module package needs to be in the same directory as MicroWorldOmics.exe.

6. Open and launch “MicroWorldOmics” in your start menu:

